

## FoxO4 (Phospho-Ser197) Polyclonal Antibody

Catalog No: #12200



Package Size: #12200-1 50ul #12200-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FoxO4 (Phospho-Ser197) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-FoxO4 (S197) Polyclonal Antibody detects endogenous levels of FoxO4 protein only when phosphorylated at S197.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human FoxO4 around the phosphorylation site of S197.
Target Name	FoxO4
Modification	Phospho
Other Names	FOXO4; AFX; AFX1; MLLT7; Forkhead box protein O4; Fork head domain transcription factor AFX1
Accession No.	Swiss-Prot: P98177NCBI Gene ID: 4303
Uniprot	P98177
GeneID	4303;
Target Species	human
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/20000

Not yet tested in other applications

## Images



Western Blot analysis of COLO205 cells using Phospho-FoxO4 (S197) Polyclonal Antibody

---

Note: This product is for in vitro research use only