

GRF-1 (Phospho-Tyr1087) Polyclonal Antibody

Catalog No: #12201

Package Size: #12201-1 50ul #12201-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GRF-1 (Phospho-Tyr1087) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-GRF-1 (Y1087) Polyclonal Antibody detects endogenous levels of GRF-1 protein only when phosphorylated at Y1087.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human GRF-1 around the phosphorylation site of Y1087.
Target Name	GRF-1
Modification	Phospho
Other Names	ARHGAP35; GRF1; GRLF1; KIAA1722; Rho GTPase-activating protein 35; Glucocorticoid receptor DNA-binding factor 1; Glucocorticoid receptor repression factor 1; GRF-1; Rho GAP p190A; p190-A
Accession No.	Swiss-Prot: Q9NRY4NCBI Gene ID: 2909
Uniprot	Q9NRY4
GeneID	2909;
Target Species	human
SDS-PAGE MW	190kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000

ELISA: 1/10000

Not yet tested in other applications

Images



Western Blot analysis of 293T cells using Phospho-GRF-1 (Y1087) Polyclonal Antibody

Note: This product is for in vitro research use only