

MEK-3 (Phospho-Ser218) Polyclonal Antibody

Catalog No: #12206

Package Size: #12206-1 50ul #12206-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEK-3 (Phospho-Ser218) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-MEK-3 (S189) Polyclonal Antibody detects endogenous levels of MEK-3 protein only when phosphorylated at S189.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MEK-3 around the phosphorylation site of S189.
Target Name	MEK-3
Modification	Phospho
Other Names	MAP2K3; MEK3; MKK3; PRKMK3; SKK2; Dual specificity mitogen-activated protein kinase kinase 3; MAP kinase kinase 3; MAPKK 3; MAPK/ERK kinase 3; MEK 3; Stress-activated protein kinase kinase 2; SAPK kinase 2; SAPKK-2; SAPKK2
Accession No.	Swiss-Prot: P46734NCBI Gene ID: 5606
Uniprot	P46734
GeneID	5606;
Target Species	human
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/5000

Not yet tested in other applications

Images

Western Blot analysis of HeLa cells using Phospho-MEK-3 (S218) Polyclonal Antibody



Note: This product is for in vitro research use only