

MEK-4 (Phospho-Thr261) Polyclonal Antibody

Catalog No: #12208

Package Size: #12208-1 50ul #12208-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEK-4 (Phospho-Thr261) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-MEK-4 (T261) Polyclonal Antibody detects endogenous levels of MEK-4 protein only when phosphorylated at T261.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MEK-4 around the phosphorylation site of T261.
Target Name	MEK-4
Modification	Phospho
Other Names	MAP2K4; JNKK1; MEK4; MKK4; PRKMK4; SEK1; SERK1; SKK1; Dual specificity mitogen-activated protein kinase kinase 4; MAP kinase kinase 4; MAPKK 4; JNK-activating kinase 1; MAPK/ERK kinase 4; MEK 4; SAPK/ERK kinase 1; SEK1; Stress-activated pro
Accession No.	Swiss-Prot: P45985NCBI Gene ID: 6416
Uniprot	P45985
GeneID	6416;
Target Species	human
SDS-PAGE MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/10000

Not yet tested in other applications

Images

Western Blot analysis of 293 cells using Phospho-MEK-4 (T261) Polyclonal Antibody



Note: This product is for in vitro research use only