

p53 (Phospho-Ser9) Polyclonal Antibody

Catalog No: #12214

Package Size: #12214-1 50ul #12214-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p53 (Phospho-Ser9) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Mk
Specificity	Phospho-p53 (S9) Polyclonal Antibody detects endogenous levels of p53 protein only when phosphorylated at S9.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human p53 around the phosphorylation site of S9.
Target Name	p53
Modification	Phospho
Other Names	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
Accession No.	Swiss-Prot: P04637NCBI Gene ID: 7157
Uniprot	P04637
GeneID	7157;
Target Species	human
SDS-PAGE MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/10000

Not yet tested in other applications

Images

Western Blot analysis of COS7 cells using Phospho-p53 (S9)
Polyclonal Antibody



Note: This product is for in vitro research use only