

Csk (Phospho-Ser364) Polyclonal Antibody

Catalog No: #12225

Package Size: #12225-1 50ul #12225-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Csk (Phospho-Ser364) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Csk (S364) Polyclonal Antibody detects endogenous levels of Csk protein only when phosphorylated at S364.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Csk around the phosphorylation site of S364.
Target Name	Csk
Modification	Phospho
Other Names	CSK; Tyrosine-protein kinase CSK; C-Src kinase; Protein-tyrosine kinase CYL
Accession No.	Swiss-Prot: P41240NCBI Gene ID: 1445
Uniprot	P41240
GeneID	1445;
Target Species	human
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/10000

Not yet tested in other applications

Images

Western Blot analysis of VEC cells using Phospho-Csk (S364)
Polyclonal Antibody



Note: This product is for in vitro research use only