

## Recombinant Human GM-CSF,CSF2 (E. coli)

Catalog No: #AP74770



Package Size: #AP74770-1 10ug #AP74770-2 50ug #AP74770-3 500ug #AP74770-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GM-CSF,CSF2 (E. coli)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04141
Uniprot	P04141
GeneID	1437;
Calculated MW	14.6kDa
Target Sequence	MAPARSPSPSTQPWEHVNAIQEARLLNLSRDAAEMNETVEVISEMFDLQEPTCLQTRLELYKQGLRGS LTK LKGPLTMMASHYKQHCPPTPETSCATQIITFESFKENLKDFLLVIPFDCWEPVQE
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Mannitol, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only