

Recombinant Human Albumin,HSA Fusion Interferon α-2b



Catalog No: #AP74771

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74771-1 10ug #AP74771-2 50ug #AP74771-3 500ug #AP74771-4 1mg

Description

Product Name	Recombinant Human Albumin,HSA Fusion Interferon α-2b
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01563
Uniprot	P01563
GeneID	3440;
Calculated MW	85.7kDa
Target Sequence	<p>DAHKSEVAHRFKDLGEEFNFKALVLIQFAQYLQQCPFEDHVKLVNEVTEFAKTCVADESAENCDKSLHTLFGDK LCTVATLRETYGEMADCCAKQEPERNECFQHKDDNPNLPRLVPEVDVMCTAFHDNEETFLKLYEYIARR HPYFYAPELLFFAKRYKAAFTCCQAADKAAACLLPKLDELDEGKASSAKQRLKQKASLQKFGERAFKAWAVA RLSQRFPKAEFAEVSKLVTDLTKVHTECCHGDLLECADRADLAKYICENQDSISSKLEKCEKPLLEKSHCIA EVENDEMPADLPSLAADFVESKDVCKNYAEAKDVFLGMFLYFYARRHPDYSVLLLRLLAKTYETTLEKCCAAA DPHECYAKVFDEFKPLVEEPQNLKQNCLEFQELGEYKFNALLVRYTKKVPQVSTPTLVEVSRNLGKVGSKC CKHPEAKRMPCAEDYLSVVLNQLCVLHEKTPVSDRVTKCCTESLVNRRPCFSALEVDETYVPKEFNAETFTFH ADICTLSEKERQIKKQATLVELVKHKPKATKEQLKAVMDDFAAFVEKCKADDKETCFAEEGKLVAAASQAAL GLCDLPQTHSLGSRRTLMLLAQMRRISLFSCLKDRHDFGFPQEEFGNQFQKAETIPVLHEMIQQIFNLFSTKDS SAAWDETLLDKFYTELYQQLNDLEACVIQGVGTETPLMKEDSILAVRKYFQRITLYLKEKKYSPCAWEVRAE IMRSFSLSTNLQESLRSKE</p>
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only