

Recombinant Human Interferon α 2B Variant (Arg46), IFN- α 2B

Catalog No: #AP74772

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74772-1 10ug #AP74772-2 50ug #AP74772-3 500ug #AP74772-4 1mg

Description

Product Name	Recombinant Human Interferon α 2B Variant (Arg46), IFN- α 2B
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01563
Uniprot	P01563
GeneID	3440;
Calculated MW	19.4kDa
Target Sequence	MCDLPQTHSLGSRRTLMLLAQMRRISLFSCLKDRHDFGFPQEEFGNQFQKAETIPVLHEMIQQIFNLFSTKDSS AAWDETLLDKFYTELYQQLNDLEACVIQGVGTETPLMKEDSILAVRKYFQRITLYLKEKKYSPCAWEVVRAEI MRSFSLSTNLQESLRSKE
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only