

## Recombinant Human Interleukin-11,IL-11

Catalog No: #AP74773



Package Size: #AP74773-1 10ug #AP74773-2 50ug #AP74773-3 500ug #AP74773-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Human Interleukin-11,IL-11   |
| Host Species    | E. coli  |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Human  |
| Accession No.   | P20809   |
| Uniprot         | P20809   |
| GeneID          | 3589;  |
| Calculated MW   | 19kDa  |
| Target Sequence | GPPPGPPRVSPDPRAELDSTVLLTRSLLDTRQLAAQLRDKFPADGDHNLDSLPTLAMSAGALGALQLPGVLT<br>RLRADLLSYLRHVQWLRRRAGGSSSLKTLEPELGTQLQARLDRLRLQLLMSRLALPQPPDPAPPLAPPSSAW<br>GGIRAAHAILGGLHLTDWAVRGLLLKTRL  |
| Formulation     | Lyophilized from a 0.2 µm filtered solution of 20mM PB, 2% Glycine, pH 7.2.  |
| Storage         | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only