

Recombinant Human Tumor Necrosis Factor  $\alpha$ , TNF $\alpha$ 

Catalog No: #AP74774



Package Size: #AP74774-1 10ug #AP74774-2 50ug #AP74774-3 500ug #AP74774-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Tumor Necrosis Factor $\alpha$ , TNF $\alpha$
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01375
Uniprot	P01375
GeneID	7124;
Calculated MW	17.3kDa
Target Sequence	MVRSSSRTPSDKPVAVHVVANPQAEGLQWLNRANALLANGVELRDNQLVVPSEGLYLIYSQVLFKGGGCP STHVLLTHTISRIAVSYQTKVNLLSAIKSPCQRETPEGAEAKPWYEPIYLGGVFQLEKGDRLSAEINRPDYLDFA ESGQVYFGIIAL
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only