

Recombinant Human Parathyroid Hormone, Parathormone, PTH (1-84)

Catalog No: #AP74775

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74775-1 10ug #AP74775-2 50ug #AP74775-3 500ug #AP74775-4 1mg

Description

Product Name	Recombinant Human Parathyroid Hormone, Parathormone, PTH (1-84)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01270
Uniprot	P01270
GeneID	5741;
Calculated MW	9kDa
Target Sequence	SVSEIQLMHNLGKHLNSMERVEWLRKKLQDVHNFVALGAPLAPRDAGSQRPRKKEDNVLVESHEKSLGEAD KADVNVLTAKSQ
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM HAc-NaAc, 150mM NaCl, 5% Mannitol, pH 4.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only