

Recombinant Human CD40 Ligand,CD40L,TNFSF5 (E. coli)



Catalog No: #AP74776

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74776-1 10ug #AP74776-2 50ug #AP74776-3 500ug #AP74776-4 1mg

Description

Product Name	Recombinant Human CD40 Ligand,CD40L,TNFSF5 (E. coli)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P29965
Uniprot	P29965
GeneID	959;
Calculated MW	16.2kDa
Target Sequence	MQKGDQNPQIAAHVISEASSKTTSVLQWAEKGYTMSNNLVTLENGKQLTVKRQGLYYIAQVTFCSNREAS SQAPFIASLCLKSPGRFRERILLRAANTHSSAKPCGQQSIHLGGVFELQPGASVFNVTDPQSQVSHGTGFTSFG LLKL
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 200mM NaCl, 0.1mM EDTA, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only