

# Recombinant Human Bone Morphogenetic Protein 2,BMP-2



Catalog No: #AP74777

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74777-1 10ug #AP74777-2 50ug #AP74777-3 500ug #AP74777-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Bone Morphogenetic Protein 2,BMP-2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P12643
Uniprot	P12643
GeneID	650;
Calculated MW	13.3kDa
Target Sequence	MQAKHKQRKRLKSSCKRHPLYVDFSDVGWNDWIVAPPGYHAFYCHGECPPFLADHLNSTNHAIVQTLVNSV NSKIPKACCVPTLSAISMLYLDENEKVVLKKNYQDMVVEGCGCR
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM HAC .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only