

Recombinant Human Interferon γ , IFN- γ (E. coli)

Catalog No: #AP74779



Package Size: #AP74779-1 10ug #AP74779-2 50ug #AP74779-3 500ug #AP74779-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human Interferon γ , IFN- γ (E. coli) |
| Host Species | E.coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | P01579 |
| Uniprot | P01579 |
| GeneID | 3458; |
| Calculated MW | 16.88kDa |
| Target Sequence | MQDPYVKEAENLKKYFNAGHSDVADNGTLFLGILKNWKEESDRKIMQSQIVSFYFKLFKNFKDDQSIQKSVETI KEDMNVKFFNSNKKKRDDFEKLTNYSVTDLNVQRKAIHELIQVMAELSPAAKTGKRKRSQMLFRGRRASQ |
| Formulation | Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.0. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only