

Recombinant Human Interleukin-15,IL-15

Catalog No: #AP74780



Package Size: #AP74780-1 10ug #AP74780-2 50ug #AP74780-3 500ug #AP74780-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-15,IL-15
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P40933
Uniprot	P40933
GeneID	3600;
Calculated MW	12.5kDa
Target Sequence	NWVNVISDLKKIEDLIQSMHIDATLYTESDVHPCKVTAMKCFLELQVISLES GDASIHDTVENLILANNSLSSN GNVTESGCKECEEELEEKNIKEFLQSFVHIVQMFINTS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only