

Recombinant Human LR3 Insulin-Like Growth Factor I,LR3-IGF-1 (MG)

Catalog No: #AP74785

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74785-1 10ug #AP74785-2 50ug #AP74785-3 500ug #AP74785-4 1mg

Description

Product Name	Recombinant Human LR3 Insulin-Like Growth Factor I,LR3-IGF-1 (MG)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05019
Uniprot	P05019
GeneID	3479;
Calculated MW	9.1kDa
Target Sequence	MFPAMPLSSLFVNGPRTLGAELVDALQFVCGDRGFYFNKPTGYGSSRRAPQTGIVDECCFRSCDLRRLEM YCAPLKPAKSA
Formulation	Lyophilized from a 0.2 µm filtered solution of 300mM HAc-NaAc, pH 6.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only