

# Recombinant Human Interferon $\alpha$ 2A Variant (Lys46),IFN- $\alpha$ 2A

Catalog No: #AP74786

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74786-1 10ug #AP74786-2 50ug #AP74786-3 500ug #AP74786-4 1mg

## Description

Product Name	Recombinant Human Interferon $\alpha$ 2A Variant (Lys46),IFN- $\alpha$ 2A
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01563
Uniprot	P01563
GeneID	3440;
Calculated MW	19.24kDa
Target Sequence	MCDLPQTHSLGSRRTLMLLAQMRKISLFSCLKDRHDFGFPQEEFGNQFQKAETIPVLHEMIQQIFNLFSTKDSS AAWDETLDDKFYTELYQQLNDLEACVIQGVGVTEPLMKEDSILAVRKYFQRITLYLKEKKYSPCAWEVVRAEI MRSFSLSTNLQESLRSKE
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only