

Recombinant Human Insulin-Like Growth Factor I, IGF-I, IGF1 (4-70)

Catalog No: #AP74789

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74789-1 10ug #AP74789-2 50ug #AP74789-3 500ug #AP74789-4 1mg

Description

Product Name	Recombinant Human Insulin-Like Growth Factor I, IGF-I, IGF1 (4-70)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05019
Uniprot	P05019
GeneID	3479;
Calculated MW	7.3kDa
Target Sequence	TLCGAELVDALQFVCGDRGFYFNKPTGYGSSSRAPQTGIVDECCFRSCDLRRLEMYCAPLKPAKSA
Formulation	Lyophilized from a 0.2 µm filtered solution of 300mM NaAc, pH 6.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only