

# Recombinant Human Insulin-Like Growth Factor I, IGF-I, IGF1



Catalog No: #AP74790

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74790-1 10ug #AP74790-2 50ug #AP74790-3 500ug #AP74790-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Insulin-Like Growth Factor I, IGF-I, IGF1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05019
Uniprot	P05019
GeneID	3479;
Calculated MW	7.6kDa
Target Sequence	GPETLCGAELVDALQFVCGDRGFYFNKPTGYGSSRRAPQTGIVDECCFRSCDLRRLRLEMYCAPLKPAKSA
Formulation	Lyophilized from a 0.2 µm filtered solution of 300mM NaAc, pH 6.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only