

Recombinant Human IL-1 Receptor Antagonist Protein,IL-1RA,IL-1F3



Catalog No: #AP74791

Orders: order@signalwayantibody.com

Package Size: #AP74791-1 10ug #AP74791-2 50ug #AP74791-3 500ug #AP74791-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL-1 Receptor Antagonist Protein,IL-1RA,IL-1F3
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P18510
Uniprot	P18510
GeneID	3557;
Calculated MW	17.26kDa
Target Sequence	MRPSGRKSSKMQAFRIWDVNQKTFYLRNNQLVAGYLQGPNVNLEEKIDVVPIEPHALFLGIHGGKMCLSCVKS GDETRLQLEAVNITDLSNRKQDKRFAFIRSDSGPTTSFESAACPGWFLCTAMEADQPVSLTNMPDEGVMVT KFYFQEDE
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM TrisHCl, 200mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only