

Recombinant Human Stem Cell Factor,SCF,c-Kit Ligand

Catalog No: #AP74792

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74792-1 10ug #AP74792-2 50ug #AP74792-3 500ug #AP74792-4 1mg

Description

Product Name	Recombinant Human Stem Cell Factor,SCF,c-Kit Ligand
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P21583
Uniprot	P21583
GeneID	4254;
Calculated MW	18.4kDa
Target Sequence	MEGICRNRVTNNVKDVTKLVANLPKDYMITLKYVPGMDVLP SHCWISEM VVQLSDSLTDLLDKFSNISEGLSNY SIIDKLVNIVDDLVECVKENS SKDLKKSFKSPEPRLFTPEEFFRIFNRSIDAFKDFVVASETSDCVV SSTLSPEKD SRVSVTKPFMLPPVA
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only