

## Recombinant Human Interleukin-8,IL-8 (Ala23-Ser99)

Catalog No: #AP74795



Package Size: #AP74795-1 10ug #AP74795-2 50ug #AP74795-3 500ug #AP74795-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Interleukin-8,IL-8 (Ala23-Ser99)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P10145
Uniprot	P10145
GeneID	3576;
Calculated MW	8.9kDa
Target Sequence	AVLPRSAKELRCQCIKTYSKPFHPKFIKELRVIESGPHCANTEIIVKLSDBGRELCLDPKENWVQRVVEKFLKRAE NS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only