

c-Fms (Phospho-Tyr809) Polyclonal Antibody

Catalog No: #12229

Package Size: #12229-1 50ul #12229-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	c-Fms (Phospho-Tyr809) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-c-Fms (Y809) Polyclonal Antibody detects endogenous levels of c-Fms protein only when phosphorylated at Y809.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human c-Fms around the phosphorylation site of Y809.
Target Name	c-Fms
Modification	Phospho
Other Names	CSF1R; FMS; Macrophage colony-stimulating factor 1 receptor; CSF-1 receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD antigen CD115
Accession No.	Swiss-Prot: P07333NCBI Gene ID: 1436
Uniprot	P07333
GeneID	1436;
Target Species	human
SDS-PAGE MW	108kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000

ELISA: 1/5000

Not yet tested in other applications

Images



Western blot analysis of 293 using p-c-Fms (Y809) antibody.

Note: This product is for in vitro research use only