

Recombinant Mouse Interleukin-1 α ,IL-1 α

Catalog No: #AP74797



Package Size: #AP74797-1 10ug #AP74797-2 50ug #AP74797-3 500ug #AP74797-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-1 α ,IL-1 α
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P01582
Uniprot	P01582
GeneID	16175;
Calculated MW	18kDa
Target Sequence	SAPYTYQSDLRYKLMKLVQRQKFMNDSLNQTIYQDVKHYLSTTWLNDLQQEVKFDMYAYSSGGDDSKYPV TLKISDSQLFVSAQGEDQPVLLKELPETPKLITGSETDLIFFWKSINSKNYFTSAAYPELFIATKEQSRVHLARGL PSMTDFQIS
Formulation	Lyophilized from a 0.2 μ m filtered solution of 50mM TrisHCl, 200mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only