

# Recombinant Mouse Fibroblast Growth Factor 1,FGF-1,FGFa (Phe16-Asp155)

Catalog No: #AP74799

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74799-1 10ug #AP74799-2 50ug #AP74799-3 500ug #AP74799-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Fibroblast Growth Factor 1,FGF-1,FGFa (Phe16-Asp155)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P61148
Uniprot	P61148
GeneID	14164;
Calculated MW	15.7kDa
Target Sequence	FNLPLGNYKKPKLLYCSNGGHFLRILPDGTVDGTRDRSDQHIQLQLSAESAGEVYIKGTETGQYLAMDTEGLL YGSQTPNEECLFLERLEENHYNTYTSKKHAEKNWVGLKKNKNGSCKRGRPRTHYGQKAILFLPLPVSSD
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 500mM NaCl, pH 6.6.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only