

Recombinant Mouse Fibroblast Growth Factor 2,FGF-2,FGFb (Met1-Ser154)

Catalog No: #AP74800

Orders: order@signalwayantibody.com

Package Size: #AP74800-1 10ug #AP74800-2 50ug #AP74800-3 500ug #AP74800-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Fibroblast Growth Factor 2,FGF-2,FGFb (Met1-Ser154)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P15655
Uniprot	P15655
GeneID	14173;
Calculated MW	17.15kDa
Target Sequence	MAASGITSLPALPEDGGAAFPPGHFKDPKRLYCKNGGFFLRIHPDGRVDGVREKSDPHVKLQLQAEERGVSISKGVCANRYLAMKEDGRLLASKCVTEECFFFERLESNNYNTYRSRKYSSWYVALKRTGQYKLGSKTGPQGQKAILFLPMSAKS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 400mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only