

Recombinant Human Fibroblast Growth Factor 2,FGF-2,FGFb (Gly132-Ser288)

Catalog No: #AP74802

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74802-1 10ug #AP74802-2 50ug #AP74802-3 500ug #AP74802-4 1mg

Description

Product Name	Recombinant Human Fibroblast Growth Factor 2,FGF-2,FGFb (Gly132-Ser288)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09038
Uniprot	P09038
GeneID	2247;
Calculated MW	16.4kDa
Target Sequence	GTMAAGSITTLPALPEDGGSGAFPPGHFKDPKRLYCKNGGFFLRIHPDGRVDGVREKSDPHIKLQLQAEERG VVSIGVCANRYLAMKEDGRLLASKCVTDECFERLESNNYNTYRSRKYTSWYVALKRTGQYKLGSKTGPG QKAILFLPMSAKS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only