

Recombinant Mouse Interleukin-22,IL-22 (E. coli)

Catalog No: #AP74803



Package Size: #AP74803-1 10ug #AP74803-2 50ug #AP74803-3 500ug #AP74803-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-22,IL-22 (E. coli)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9JJY9
Uniprot	Q9JJY9
GeneID	50929;
Calculated MW	16.7kDa
Target Sequence	MLPVNTRCKLEVSNFQQPYIVNRTFMLAKEASLADNNTDVRLIGEKLFRGVS AKDQC CYLMKQVLNFTLEDVLL PQSDRFQPYMQEVV PFLTKLSNQLSSCHISGDDQNIQKNVRRLLKETVKKLGESGEIKAIGELDLLFMSLRNACV
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only