

Recombinant Mouse Growth Hormone,GH

Catalog No: #AP74804



Package Size: #AP74804-1 10ug #AP74804-2 50ug #AP74804-3 500ug #AP74804-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Growth Hormone,GH
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P06880
Uniprot	P06880
GeneID	14599;
Calculated MW	21.9kDa
Target Sequence	MFPAMPLSSLFSNAVLRAQHLHQLAADTYKEFERAYIPEGQRYSIQNAQA AFCFSETIPAPTGKEEAQQR TDM ELLRFSLLLIQSWLGPVQFLSRIFTNSLMFGTSDRVYEKLDLEEGIQALMQELEDGSPRVGQILKQTYDKFDA NMRSDALLKNYGLLSCFKDLHKAETYL RVMKCRRFVESSCAF
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM TrisHCl, 500mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only