

Recombinant Human Myelin Oligodendrocyte Glycoprotein,MOG



Catalog No: #AP74807

Orders: order@signalwayantibody.com

Package Size: #AP74807-1 10ug #AP74807-2 50ug #AP74807-3 500ug #AP74807-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Myelin Oligodendrocyte Glycoprotein,MOG
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16653
Uniprot	Q16653
GeneID	4340;
Calculated MW	15.2kDa
Target Sequence	MGQFRVIGPRRHPICALVGDEVELPCRISPGKNATGMEVGWYRPPFSRVVHLYRNGKDDGDQAPEYRGRTE LLKDAIGEGKVTLRIRNVRFSDEGGFTCFRDHSYQEEAAMELKVEDPFYWVSPGHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM HAC-NaAc, 150mM NaCl, pH 4.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only