

# Recombinant Human C-X-C Motif Chemokine 10,CXCL10

Catalog No: #AP74809

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP74809-1 10ug #AP74809-2 50ug #AP74809-3 500ug #AP74809-4 1mg

## Description

Product Name	Recombinant Human C-X-C Motif Chemokine 10,CXCL10
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02778
Uniprot	P02778
GeneID	3627;
Calculated MW	8.6kDa
Target Sequence	VPLSRTVTRCTCISISNQPVNPRSLEKLEIIPASQFCPRVEIIATMKKKGEKRCLNPESKAIKNLLKAVSKERSKRS P
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only