

Recombinant Human C-X-C Motif Chemokine 14,CXCL14



Catalog No: #AP74810

Orders: order@signalwayantibody.com

Package Size: #AP74810-1 10ug #AP74810-2 50ug #AP74810-3 500ug #AP74810-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human C-X-C Motif Chemokine 14,CXCL14
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95715
Uniprot	O95715
GeneID	9547;
Calculated MW	9.4kDa
Target Sequence	SKCKCSRKGPKIRYSDVKKLEMKPKYPHCEEKMMVIITTKSVSRYRGQEHCLHPKLQSTKRFIKWYNWNEKRR VYEE
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 1M NaCl, pH 8.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only