

PPAR- γ (Phospho-Ser112) Polyclonal Antibody

Catalog No: #12231

Package Size: #12231-1 50ul #12231-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPAR- γ (Phospho-Ser112) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-PPAR- γ (S112) Polyclonal Antibody detects endogenous levels of PPAR- γ protein only when phosphorylated at S112.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human PPAR- γ around the phosphorylation site of S112.
Target Name	PPAR- γ
Modification	Phospho
Other Names	PPARG; NR1C3; Peroxisome proliferator-activated receptor gamma; PPAR-gamma; Nuclear receptor subfamily 1 group C member 3
Accession No.	Swiss-Prot: P37231NCBI Gene ID: 5468
Uniprot	P37231
GeneID	5468;
Target Species	human
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000


ELISA: 1/10000

Not yet tested in other applications

Images



Western Blot analysis of PC-3 cells using Phospho-PPAR- γ (S112) Polyclonal Antibody



Western Blot analysis of various cells using Phospho-PPAR- γ (S112) Polyclonal Antibody

Note: This product is for in vitro research use only