

Recombinant Human Interleukin-33,IL-33

Catalog No: #AP74831



Package Size: #AP74831-1 10ug #AP74831-2 50ug #AP74831-3 500ug #AP74831-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-33,IL-33
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95760
Uniprot	O95760
GeneID	90865;
Calculated MW	20.29kDa
Target Sequence	MSITGISPITEYLASLSTYNDQSITFALEDESYEIYVEDLKKDEKKDKVLLSYYESQHPSPNESGDGVDGKMLMVT LSPTKDFWLHANNKEHSVELHKCEKPLPDQAFFVLHNMHSNCVSECKTDPGVFIGVKNHLALIKVDSSENL CTENILFKLSET
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only