

Recombinant Human Annexin A1,ANXA1

Catalog No: #AP74832



Package Size: #AP74832-1 10ug #AP74832-2 50ug #AP74832-3 500ug #AP74832-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Annexin A1,ANXA1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04083
Uniprot	P04083
GeneID	301;
Calculated MW	38.58kDa
Target Sequence	AMVSEFLKQAWFIENEEQEYVQTVKSSKGGPGSAVSPYPTFNPSSDVAALHKAIMVKGVDIATIIDILTKRNN QRQIQIAAYLQETGKPLDETLKKALTGHLEEVVLLALLKTPAQFDADELRAAMKGLGTDEDTLIEILASRTNKEIR DINRVYREELKRD LAKDITSDTSGDFRNALLSLAKGDRSEDFGVNEDLADSDARALYEAGERRKGTDVNVFNTI LTTRSYPQLRRVFQKYTKYSKHD MNKVL DLELKGDIKCLTAIVKCATSKPAFFAEK LHQAMKGVGTRHKALIRI MVSRS EIDMNDIKAFYQKMYGISLCQAILDETKGDYEKILVALCGGN
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only