

Recombinant Human Acyl-Protein Thioesterase 1,APT-1 (N-6His)

Catalog No: #AP74835

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74835-1 10ug #AP74835-2 50ug #AP74835-3 500ug #AP74835-4 1mg

Description

Product Name	Recombinant Human Acyl-Protein Thioesterase 1,APT-1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75608
Uniprot	O75608
GeneID	10434;
Calculated MW	26.83kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMCGNNMSTPLPAIVPAARKATAAVIFLHGLGDTGHGWAEAFAGIRSSHIKYIC PHAPVRPVTLNMMNVAMPSWFDIIGLSPDSQEDESIGKQAENIKALIDQEVKNGIPSNRIILGGFSQGGALSLYT ALTTQQKLAGVTALSCWLPLRASFPQGPiGGANRDISILQCHGDCDPLVPLMFGSLTVEKLTLPANVTFKT YEGMMHSSCQEMMDVKQFIDKLLPPID
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only