

Recombinant Human Ciliary Neurotrophic Factor, CNTF



Catalog No: #AP74836

Orders: order@signalwayantibody.com

Package Size: #AP74836-1 10ug #AP74836-2 50ug #AP74836-3 500ug #AP74836-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ciliary Neurotrophic Factor, CNTF
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P26441
Uniprot	P26441
GeneID	1270;
Calculated MW	22.93kDa
Target Sequence	AFTEHSPLTPHRRDLCRSIWLARKIRSDLTALTESYVKHQGLNKNINLDSADGMPVASTDQWSELTEAERLQ ENLQAYRTFHVLLARLLEDQQVHFTPTGDFHQAIHTLLLQVAAFAYQIEELMILLEYKIPRNEADGMPINVGDG GLFEKKLWGLKVLQELSQWTVRSIHDLRFISSHQTGIPARGSHYIANNKKM
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM Cysteine, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only