

Recombinant Human Apurinic-Apyrimidinic Endonuclease 1,APE

Catalog No: #AP74839

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74839-1 10ug #AP74839-2 50ug #AP74839-3 500ug #AP74839-4 1mg

Description

Product Name	Recombinant Human Apurinic-Apyrimidinic Endonuclease 1,APE
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P27695
Uniprot	P27695
GeneID	328;
Calculated MW	35.62kDa
Target Sequence	MGPKRGKKGAVAEDGDELRTPEAKKSKTAAKKNCKEAAAGEGPALYEDPPDQKTSPSGKPATLKICSWNVD GLRAWIKKKGLDWVKEEAPDILCLQETKCSENKLPALQELPGLSHQYWSAPSDKEGYSGVGLLSRQCPLKV SYGIGEEEHQDEGRVIVAEFDSFVLVTAYVPNAGRGLVRLEYRQRWDEAFRKFLLKGLASRKPLVLCGDLNVA HEEIDLRNPKGNKKNAGFTPQERQGFCELLQAVPLADSFRLHPNTPYAYTFWTYMMNARSKNVGWRLDYF LLSHSLLPALCDSKIRSKALGSDHCPITLYLAL
Formulation	Supplied as a 0.2 µm filtered solution of 10mM HEPES, 100mM KCl, 50% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only