

Recombinant Human Bcl-2 Related Protein A1,BCL2A1 (C-6His)

Catalog No: #AP74841

Orders: order@signalwayantibody.com

Package Size: #AP74841-1 10ug #AP74841-2 50ug #AP74841-3 500ug #AP74841-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Bcl-2 Related Protein A1,BCL2A1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16548
Uniprot	Q16548
GeneID	597;
Calculated MW	18.52kDa
Target Sequence	MTDCEFGYIYRLAQDYLCVQLQIPQPGSGPSKTSRVLQNVAFSVQKEVEKNLKSCLDNVNVVSVDTARTLFNQ VMEKEFEDGIINWGRIVTIFAFEGILIKLLRQQIAPDVDTYKEISYFVAEFIMNNTGEWIRQNGGWENGWVKKFE PKSLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 1mM DTT, 10% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only