

Recombinant Human Bcl-2-Like Protein 1, BCL2L1 (C-6His)

Catalog No: #AP74842

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74842-1 10ug #AP74842-2 50ug #AP74842-3 500ug #AP74842-4 1mg

Description

Product Name	Recombinant Human Bcl-2-Like Protein 1, BCL2L1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q07817
Uniprot	Q07817
GeneID	598;
Calculated MW	24.85kDa
Target Sequence	MSQSNRELVVDFLSYKLSQKGYSSWSQFSDVEENRTEAPEGTESEMETPSAINGNPSWHLADSPAVNGATGH SSSLDAREVIPMAAVKQALREAGDEFELRYRRAFSDLTSQLHITPGTAYQSFEQVVNELFRDGVNWGRIVAFF SFGGALCVESVDKEMQVLVSRIAAMATYLNHDHLEPWIQENGGWDTFVELYGNAAAESRKGQERFNRLEH HHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM HEPES, 50mM KCl, 20% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only