

Recombinant Human Carbonic Anhydrase 1,CA1 (C-6His)

Catalog No: #AP74845

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74845-1 10ug #AP74845-2 50ug #AP74845-3 500ug #AP74845-4 1mg

Description

Product Name	Recombinant Human Carbonic Anhydrase 1,CA1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00915
Uniprot	P00915
GeneID	759;
Calculated MW	29.93kDa
Target Sequence	ASPDWGYDDKNGPEQWSKLYPIANGNNQSPVDIKTSETKHDTSLKPISVSYNPATAKEIINVGHFSFHVNFEDN DNRSVLKGGPFSDSYRLFQFHFHWGSTNEHGSEHTVDGVKYSaelhVAHWNSAKYSSLAEAASKADGLAVI GVLmkvGEANPKLQKVLdalQAiKTKGRAPFTNFDpstllPSSldFWtYPgSLThPPlyESVtWiICKESISVS SEqlAQFRSllSNveGDNAVPMQHNNRPTQPLKGRtVRASfLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only