

Recombinant Human Carbonic Anhydrase 13,CA13 (C-6His)

Catalog No: #AP74846

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74846-1 10ug #AP74846-2 50ug #AP74846-3 500ug #AP74846-4 1mg

Description

Product Name	Recombinant Human Carbonic Anhydrase 13,CA13 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N1Q1
Uniprot	Q8N1Q1
GeneID	377677;
Calculated MW	30.51kDa
Target Sequence	MSRLSWGYPREHNGPIHWKEFFPIADGDQQSPIEIKTKEVKYDSSLRPLSIKYPSSAKIISNSGHSFNVDFFDT ENKSVLRGGPLTGSYRLRQVHLHWGSADDHGSEHIVDGVSYAAELHVVHWNSDKYPSFVEAAHEPDGLAVL GVFLQIGEPNSQLQKITDTLDSIKEKGKQTRFTNFDLLSLLPPSWDYWTYPGSLTVPPLLESVTWIVLKQPINISS QQLAKFRSLLCTAEGEAAAFLVSNHRPPQPLKGRKVRASFHLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only