

Recombinant Human Carbonic Anhydrase 4,CA4 (C-6His)

Catalog No: #AP74847

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74847-1 10ug #AP74847-2 50ug #AP74847-3 500ug #AP74847-4 1mg

Description

Product Name	Recombinant Human Carbonic Anhydrase 4,CA4 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P22748
Uniprot	P22748
GeneID	762;
Calculated MW	31.43kDa
Target Sequence	MAESHWCYEVQAESSNYPCLVPVKWGGNCQKDRQSPINIVTTKAKVDKKLGRFFFSGYDKKQTWTVQNNGH SVMMLLENKASISGGGLPAPYQAKQLHLHWSDLPLYKGSEHSLDGEHFAMEMHIVHEKEKGTSRNVKEAQDP EDEIAVLAFLEAGTQVNEGFQPLVEALSNIPKPEMSTTMAESSLLDLLPKEEKLRYFRYLGSLLTPTCDEKV VWTVFREPIQLHREQILAFSQKLYDKEQTVSMKDNVRPLQQLGQRTVIKLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only