

Recombinant Human Carbonic Anhydrase 8,CA8 (C-6His)



Catalog No: #AP74850

Orders: order@signalwayantibody.com

Package Size: #AP74850-1 10ug #AP74850-2 50ug #AP74850-3 500ug #AP74850-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Carbonic Anhydrase 8,CA8 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P35219
Uniprot	P35219
GeneID	767;
Calculated MW	34.04kDa
Target Sequence	ADLSFIEDTVAFPEKEEDEEEEEEGVEWGYEEGVEWGLVFPDANGEYQSPINLNSREARYDPSLLDVRLSPN YVVCRDCEVTNDGHTIQVILKSKSVLSGGPLPQGHEFELYEVRFHWGRENQRGSEHTVNFKAFFMELHLIHW NSTLFGSIDEAVGKPHGIAIIALFVQIGKEHVGLKAVTEILQDIQYK GKSKTIPCFNPNTLLPDLLRDYWVYEGSL TIPPCSEGVTWILFRYPLTISQLQIEEFRRLRTHVKGAEELVEGCDGILGDNFRPTQPLSDRVIRAAAFQLEHHHHH H
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 500mM NaCl, 1mM DTT, pH 8.5 .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only