

# Recombinant Human C-C Motif Chemokine 26,CCL26 (24-94)



Catalog No: #AP74851

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP74851-1 10ug #AP74851-2 50ug #AP74851-3 500ug #AP74851-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human C-C Motif Chemokine 26,CCL26 (24-94)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y258
Uniprot	Q9Y258
GeneID	10344;
Calculated MW	8.53kDa
Target Sequence	MTRGSDISKTCFCFYSHKPLPWTWVRSYEFTSNCSQRAVIFTTKRGKKVCTHPRKKWVQKYISLLKTPKQL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 1mM EDTA, 20% Glycerol, pH 9.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only