

Recombinant Human E3 Ubiquitin-Protein Ligase CHIP,CHIP



Catalog No: #AP74853

Orders: order@signalwayantibody.com

Package Size: #AP74853-1 10ug #AP74853-2 50ug #AP74853-3 500ug #AP74853-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human E3 Ubiquitin-Protein Ligase CHIP,CHIP
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UNE7
Uniprot	Q9UNE7
GeneID	10273;
Calculated MW	34.86kDa
Target Sequence	MKGKKEKEGGARLGAGGGSPKSPSAQELKEQGNRLFVGRKYPEAAACYGRAITRNPLVAVYYTNRALCYL KMQQHEQALADCRRALELDGQSVKAHFFLGQCQLEMESYDEAIANLQRAYSLAKEQRLNFGDDIPSALRIAKK KRWNSIEERRIHQESLHSYLSRLIAAERERELEECQRNHEGDEDDSHVRAQQACIEAKHDKYMADMDELFS QVDEKRKKRDIPDYLCGKISFELMREPCITPSGITYDRKDIEEHLQRVGHFDPVTRSPLTQEQQLIPNLAMKEVID AFISENGWVEDY
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only