

## Recombinant Human Cornulin (N-6His)

Catalog No: #AP74855



Package Size: #AP74855-1 10ug #AP74855-2 50ug #AP74855-3 500ug #AP74855-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Cornulin (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UBG3
Uniprot	Q9UBG3
GeneID	49860;
Calculated MW	17.45kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMPQLLQNINGIIEAFRRYARTEGNCALTRGELKRLLEQEFADVIVKPHDPAT VDEVLRLLDEDHTGTVEFKEFLVLFKVAQACFKTLESAEGACGSQESGSLHSGASQELGEGQQRSGTEVGR AGKGQHYEGSSHRQS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only