

Recombinant Human C-X-C Motif Chemokine 12,CXCL12,SDF-1 (22-89)

Catalog No: #AP74858

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74858-1 10ug #AP74858-2 50ug #AP74858-3 500ug #AP74858-4 1mg

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human C-X-C Motif Chemokine 12,CXCL12,SDF-1 (22-89) |
| Host Species | E.coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | P48061 |
| Uniprot | P48061 |
| GeneID | 6387; |
| Calculated MW | 8.12kDa |
| Target Sequence | KPVLSYRCPGRFFESHVARANVKHLKILNTPNCALQIVARLKNNNRQVCIDPKLKWIQEYLEKALNK |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only