

Recombinant Human β -Defensin 1,DEFB1

Catalog No: #AP74861



Package Size: #AP74861-1 10ug #AP74861-2 50ug #AP74861-3 500ug #AP74861-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human β -Defensin 1,DEFB1 |
| Host Species | E.coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | P60022 |
| Uniprot | P60022 |
| GeneID | 1672; |
| Calculated MW | 5.07kDa |
| Target Sequence | GNFLTGLGHRSDHYNCVSSGGQCLYSACPIFTKIQGTCTYRGKAKCCK |
| Formulation | Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 130mM NaCl, pH 7.4. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only